

Sterile Polynit Heatseal Low Endotoxin Wipes

The presence of pyrogens is a critical safety concern for injectable products and medical devices as products contaminated with pyrogens can pose a life-threatening risk to patients. Unlike viable microbial contaminants which can be destroyed by various sterilization techniques, pyrogens are difficult to remove and deactivate. Bacterial endotoxin, specifically from the outer membranes of gram negative bacteria, are the most common pharmaceutical pyrogens.

Using low endotoxin wipes during manufacturing reduces the potential for pyrogen contamination in the most critical of applications.

Sterile Polynit Heatseal Low Endotoxin wipes are manufactured from high quality 100% knitted polyester fabric with sealed edges and produce low levels of particles and fibers. The wipes have a guaranteed endotoxin level of less than 1 EU per wipe and are validated sterile to 10⁻⁶ SAL, making them ideal for use in product contact areas.

The wipes are half folded and triple packaged in linear tear outer bags for ease of transfer into ISO Class 5 cleanrooms. This smaller packaging takes up less space making it ideal for use in isolators and RABS. The small quantity of wipes per package can be used during one cleaning session, eliminating waste.

Features	Benefits
100% laundered knitted polyester	Very low levels of particles and fibersGood abrasion and chemical resistance
Laser sealed edges	Reduces fiber generation from the edge of the wipe
Guaranteed endotoxin levels of <1 EU per wipe	Suitable for use in product contact areas
Validated sterile or gamma irradiated to >25kGY	• Suitable for use in ISO Class 5 (Grades A/B) cleanrooms
Small number of wipes per bag	• Eliminates waste, all wipes will be used in one session
Triple packed in linear tear bags	Easy to open and transfer into isolator and RABS

Part No.	Description	Packaging
LWLE0001	Sterile Polynit Heatseal Low Endotoxin Wipes, 9" x 9" (23x23cm), Half folded, triple bagged	10/bag; 50 bags/case
LWLE0002	Sterile Polynit Heatseal Low Endotoxin Wipes, 12" x 12" (30.5x30.5cm), Half folded, triple bagged	10/bag; 36 bags/case

PDSW014 070921 Ref: 9649 Copyright © 2021 Contec, Inc. All rights reserved.







Product Information

Material	100% polyester	
Construction	Interlock knit, no-run	
Packaging Materials	Outer bags (OB1, OB2), low density polyethylene (LDPE) Case (CS), corrugated fiberboard (PAP)	
Environment	ISO 3-8 Grade A/B	



Technical Data

Attribute (units)	Typical Value	Test Method
Basis weight; nominal (g/m²)	140	Contec Method
Sorbent capacity; (mL/m²)	503	IEST-RP-CC004.3, Sec. 8.1
Sorptive rate; (seconds)	<1	
Non-volatile residue, NVR		IEST-RP-CC004.3, Sec. 7.1.1
In deionized water; (g/m²)	0.005	
In isopropanol; (g/m²)	0.008	
Specific ions		IEST-RP-CC004.3, Sec. 7.2.2
Sodium; (ppm)	0.040	
Chloride; (ppm)	0.048	
Particles, readily releasable		
$P \ge 0.5 \mu m; (x10^6/m^2)$	6.17	IEST-RP-CC004.3, Sec. 6.1.3
Fibers > $100 \mu m$; (x10 ³ /m ²)	0.113	IEST-RP-CC004.3, Sec. 6.2.2



Packaging

	EA/OB1	OB1/OB2	OB2/OB3	OB3/CS	EA/CS
LWLE0001	10	1	1	50	500
LWLE0002	10	1	1	36	360

EA = Each, OB1 = Outer Bag 1, OB2 = Outer Bag 2, OB3 = Outer bag 3, CS = Case

Recycle Symbol Key

	, ,
PET	
HDPE	22.5
LDPE	4.5
PP	
PAP	20

Notes

a) The data shown are typical values and should not be used as product specifications.
b) Valid product comparisons may only be obtained through side-by-side testing in the same test facility, under similar conditions.

c) Current and/or comparison data may be available. Please contact a Contec sales representative for details.

Test Methods:

1) CTM Contec Test Method

2) IEST-RP-CC004.3 Evaluating Wiping Materials Used in Cleanroom and Other Controlled Environments, Institute of environmental Sciences and Technology, Rolling Meadows IL.