

SUITE Touch-less products:

- Eliminate operator glove handling of both clean pads and soiled pads
- Provides consistent and liberal application of disinfectant
- Incorporates Suite engineered fiber resulting in more square feet of application per prepped pad
- Significantly reduces time spent changing pads
- Are non-particulating and snag resistant



Product: SUITE Touch-less PD

Description: Absorbent pad mop with nylon loop back, 17" x 5". For applying disinfectants on walls, ceilings and floors in cleanrooms.

Part #: Y-TL-517-25

Packaging: 5 mops per pack, 25 mops per case.

Handle:

Part #: Y-9966-M

Description: Ergonomic autoclavable anodized aluminum and plastic extendable handle.1/bag

Frame:

Part#: Y-10516-PLH

Description: Anodized Aluminum autoclavable 16" Suite TL frame with plastic lockable hinge.1/bag

Strips:

Part #: Y-10516-50

Description: Nylon hook strip replacements 2/bag, 10/pk, 5pk/cs



Testing:

<u>Test</u>	<u>Method</u>	<u>Result</u>
<u>Particles</u>	Helmke Drum	7.8 particles @ 0.5 um
<u>Metals and Ions</u>	IEST-RP-CC004.3	0.959
<u>FTIR for Silicone</u>	IEST-RP-CC004.3	None detected
<u>Sorbency Capacity</u>	IEST-RP-CC004.3	793 mL/m ²

Reduction and Transfer: Wiping with Suite material has been demonstrated to remove an average of 98%+ of surface bacterial endotoxin. In addition, the carry-over (cross contamination potential) from used Suite-Material is demonstrated to be less than 0.02%.

*Above testing performed on a 12x12 piece of Suite material

Rev: 14-01